




2.1.2 Steril-Gewindestutzen STERI-Connect®

aseptic threaded connecting piece (male part)



zusätzliche Oberflächenqualitäten analog Seite | additional surface qualities analog page 7

weitere Abmessungen auf Anfrage | further dimensions on request

DIN	Abmessung dimension	A	B	C	E (M)	MB 1.4435			
DN006	8.0 x 1.0	8.0	6.0	36.0	16 x 1.5	705205			
DN008	10.0 x 1.0	10.0	8.0	36.0	18 x 1.5	704503			
DN010	12.0 x 1.0	12.0	10.0	37.0	22 x 1.5	700059			
DN010	13.0 x 1.5	13.0	10.0	37.0	22 x 1.5	709523			
DN015	19.0 x 1.5	19.0	16.0	37.0	30 x 1.5	703799			
DN020	23.0 x 1.5	23.0	20.0	39.0	36 x 2.0	700064			
DN025	29.0 x 1.5	29.0	26.0	44.0	42 x 2.0	700066			
DN032	35.0 x 1.5	35.0	32.0	45.0	52 x 2.0	700068			
DN040	41.0 x 1.5	41.0	38.0	46.0	56 x 2.0	700070			
DN050	53.0 x 1.5	53.0	50.0	47.0	68 x 2.0	700072			
DN065	70.0 x 2.0	70.0	66.0	64.0	90 x 3.0	700075			
DN080	85.0 x 2.0	85.0	81.0	64.0	100 x 3.0	700078			
DN100	104.0 x 2.0	104.0	100.0	66.0	130 x 4.0	700061			

ISO	Abmessung dimension	A	B	C	E (M)	MB 1.4435			
DN008	13.5 x 1.6	13.5	10.3	37.0	22 x 1.5	700077			
DN010	17.2 x 1.6	17.2	14.0	37.0	26 x 1.5	700060			
DN015	21.3 x 1.6	21.3	18.1	37.0	30 x 1.5	700063			
DN020	26.9 x 1.6	26.9	23.7	39.0	36 x 2.0	700065			
DN025	33.7 x 2.0	33.7	29.7	44.0	42 x 2.0	700067			
DN032	42.4 x 2.0	42.4	38.4	45.0	52 x 2.0	700069			
DN040	48.3 x 2.0	48.3	44.3	46.0	56 x 2.0	700071			
DN050	60.3 x 2.0	60.3	56.3	47.0	68 x 2.0	700073			
DN065	76.1 x 2.0	76.1	72.1	64.0	90 x 3.0	700076			
DN080	88.9 x 2.3	88.9	84.3	64.0	100 x 3.0	700079			
DN100	114.3 x 2.3	114.3	109.7	66.0	130 x 4.0	700062			