




2.1.3 Steril-Bundstutzen STERI-Connect®

aseptic liner



weitere Abmessungen auf Anfrage | further dimensions on request

zusätzliche Oberflächenqualitäten analog Seite | additional surface qualities analog page 7

DIN	Abmessung dimension	A	B	C	D	MB 1.4435			
DN006	8.0 x 1.0	8.0	6.0	33.0	13.0	705204			
DN008	10.0 x 1.0	10.0	8.0	33.0	15.0	700033			
DN010	12.0 x 1.0	12.0	10.0	33.0	19.0	700102			
DN010	13.0 x 1.5	13.0	10.0	33.0	19.0	709767			
DN015	19.0 x 1.5	19.0	16.0	33.0	27.0	700104			
DN020	23.0 x 1.5	23.0	20.0	33.0	33.0	700105			
DN025	29.0 x 1.5	29.0	26.0	39.0	39.0	700106			
DN032	35.0 x 1.5	35.0	32.0	39.0	49.0	700107			
DN040	41.0 x 1.5	41.0	38.0	39.0	53.0	700108			
DN050	53.0 x 1.5	53.0	50.0	39.0	65.0	700109			
DN065	70.0 x 2.0	70.0	66.0	50.0	85.0	700110			
DN080	85.0 x 2.0	85.0	81.0	50.0	95.0	700111			
DN100	104.0 x 2.0	104.0	100.0	50.0	124.0	700103			

ISO	Abmessung dimension	A	B	C	D	MB 1.4435			
DN008	13.5 x 1.6	13.5	10.3	33.0	19.0	700035			
DN010	17.2 x 1.6	17.2	14.0	33.0	23.0	700025			
DN015	21.3 x 1.6	21.3	18.1	33.0	27.0	700027			
DN020	26.9 x 1.6	26.9	23.7	33.0	33.0	700028			
DN025	33.7 x 2.0	33.7	29.7	39.0	39.0	700029			
DN032	42.4 x 2.0	42.4	38.4	39.0	49.0	700030			
DN040	48.3 x 2.0	48.3	44.3	39.0	53.0	702499			
DN050	60.3 x 2.0	60.3	56.3	39.0	65.0	700032			
DN065	76.1 x 2.0	76.1	72.1	50.0	85.0	700034			
DN080	88.9 x 2.3	88.9	84.3	50.0	95.0	700036			
DN100	114.3 x 2.3	114.3	109.7	50.0	124.0	700026			