



**Zulassungen / Konformitäten**

**EPDM (peroxidisch vernetzt)**

FDA konform (§ 21 CFR 177.2600)  
USP Class VI  
3 A Sanitary Standard, Number 18-03 (Class II)  
DVGW (W270 und W534)  
KTW  
EN 681-1 WB

**FKM/FEP ummantelt**

FDA konform (§ 21 CFR 177.1550)  
USP Class VI  
3 A Sanitary Standard, Number 18-03 (Class I)

**weitere Werkstoffe auf Anfrage**

**Verpackung**

einzeln verpackt

**Approvals / Conformities**

**EPDM (crosslinked peroxidic)**

FDA compliant (§ 21 CFR 177.2600)  
USP Class VI  
3 A Sanitary Standard, Number 18-03 (Class II)  
DVGW (W270 and W534)  
KTW  
EN 681-1 WB

**FKM/FEP encapsulated**

FDA compliant (§ 21 CFR 177.1550)  
USP Class VI  
3 A Sanitary Standard, Number 18-03 (Class I)

**further materials on request**

**Packaging**

wrapped separately



weitere Abmessungen auf Anfrage | further dimensions on request

DIN	Abmessung dimension						EPDM	FKM/ FEP *		
DN006	8.0 x 1.0						705218	XXX		
DN008	10.0 x 1.0						700970	XXX		
DN010	12.0 x 1.0						700954	705221		
DN010	13.0 x 1.5						700954	705221		
DN015	19.0 x 1.5						700958	705222		
DN020	23.0 x 1.5						700960	705223		
DN025	29.0 x 1.5						700962	705224		
DN032	35.0 x 1.5						700964	705225		
DN040	41.0 x 1.5						700966	705226		
DN050	53.0 x 1.5						700968	705227		
DN065	70.0 x 2.0						700971	705228		
DN080	85.0 x 2.0						700975	705229		
DN100	104.0 x 2.0						700956	704472		

ISO	Abmessung dimension						EPDM	FKM/ FEP *		
DN008	13.5 x 1.6						700974	705230		
DN010	17.2 x 1.6						700955	705231		
DN015	21.3 x 1.6						700959	705232		
DN020	26.9 x 1.6						700961	705233		
DN025	33.7 x 2.0						700963	705234		
DN032	42.4 x 2.0						700965	705235		
DN040	48.3 x 2.0						700967	705236		
DN050	60.3 x 2.0						700969	705237		
DN065	76.1 x 2.0						700972	705238		
DN080	88.9 x 2.3						700976	705239		
DN100	114.3 x 2.3						700957	705240		

\*) ummantelt | encapsulated