


hexagonal screw set | Sechskantschraubensatz 2.8



Legende

C Gewindelänge
D Gewinde

key

C length of thread
D thread

Sechskantschraubensatz für Flanschverbindung DIN 11864-2-A-N hexagonal screw set for aseptic flange connection DIN 11864-2-A-N

Abmessung dimension DIN / ISO / ASME	C Schraube	D screw	D Mutter nut	Anzahl number	Art.-Nr. item No
010-050 / 008-050 / ½"-2"	30	M8	M8	4	700944
065 / --- / 2½"	30	M8	M8	8	700945
--- / 065 / 3"	30	M10	M10	8	700947
080 / 080 / 4"	35	M10	M10	8	700948
100 / 100 / ---	40	M10	M10	8	700949

Sechskantschraubensatz für Flanschverbindung STERI-Connect® hexagonal screw set for aseptic flange connection STERI-Connect®

Abmessung dimension DIN / ISO	C Schraube	D screw	D Mutter nut	Anzahl number	Art.-Nr. item No
008-040 / 008-040	30	M8	M8	4	700944
050 / 050	35	M8	M8	4	700946
065-080 / 065-080	40	M10	M10	4	700949
100 / 100	45	M10	M10	4	700950

Normen | standards

Schraube screw	DIN 933				
Mutter nut	DIN 934				