




2.5.4 Steril-Blindclamps -VS- STERI-Connect®

aseptic blind clamp -VS-



zusätzliche Oberflächenqualitäten analog Seite | additional surface qualities analog page 42

weitere Abmessungen auf Anfrage | further dimensions on request

| DIN | Abmessung dimension | C | E | | | | MB 1.4435 | | | |
|-------|---------------------|------|---------|--|--|--|-----------|--|--|--|
| DN006 | 8.0 x 1.0 | 7.5 | Ø 25.0 | | | | 700154 | | | |
| DN008 | 10.0 x 1.0 | 7.5 | Ø 25.0 | | | | 700156 | | | |
| DN010 | 12.0 x 1.0 | 7.5 | Ø 25.0 | | | | 700146 | | | |
| DN010 | 13.0 x 1.5 | 7.5 | Ø 25.0 | | | | 709840 | | | |
| DN015 | 19.0 x 1.5 | 7.5 | Ø 25.0 | | | | 700148 | | | |
| DN020 | 23.0 x 1.5 | 9.0 | Ø 50.5 | | | | 700149 | | | |
| DN025 | 29.0 x 1.5 | 9.0 | Ø 50.5 | | | | 700150 | | | |
| DN032 | 35.0 x 1.5 | 9.0 | Ø 50.5 | | | | 700151 | | | |
| DN040 | 41.0 x 1.5 | 9.0 | Ø 64.0 | | | | 700152 | | | |
| DN050 | 53.0 x 1.5 | 9.0 | Ø 77.5 | | | | 700153 | | | |
| DN065 | 70.0 x 2.0 | 9.0 | Ø 91.0 | | | | 700155 | | | |
| DN080 | 85.0 x 2.0 | 11.0 | Ø 106.0 | | | | 700157 | | | |
| DN100 | 104.0 x 2.0 | 11.0 | Ø 119.0 | | | | 700147 | | | |

| ISO | Abmessung dimension | C | E | | | | MB 1.4435 | | | |
|-------|---------------------|------|---------|--|--|--|-----------|--|--|--|
| DN008 | 13.5 x 1.6 | 7.5 | Ø 25.0 | | | | 704971 | | | |
| DN010 | 17.2 x 1.6 | 7.5 | Ø 25.0 | | | | 705281 | | | |
| DN015 | 21.3 x 1.6 | 7.5 | Ø 50.5 | | | | 705282 | | | |
| DN020 | 26.9 x 1.6 | 9.0 | Ø 50.5 | | | | 705283 | | | |
| DN025 | 33.7 x 2.0 | 9.0 | Ø 50.5 | | | | 705284 | | | |
| DN032 | 42.4 x 2.0 | 9.0 | Ø 50.5 | | | | 705285 | | | |
| DN040 | 48.3 x 2.0 | 9.0 | Ø 64.0 | | | | 705286 | | | |
| DN050 | 60.3 x 2.0 | 9.0 | Ø 77.5 | | | | 705287 | | | |
| DN065 | 76.1 x 2.0 | 9.0 | Ø 91.0 | | | | 705288 | | | |
| DN080 | 88.9 x 2.3 | 11.0 | Ø 106.0 | | | | 705289 | | | |
| DN100 | 114.3 x 2.3 | 11.0 | Ø 130.0 | | | | 705290 | | | |