



weitere Abmessungen auf Anfrage | further dimensions on request

DIN	Abmessung dimension	A	B	C	E	MB 1.4404
DN010	13.0 x 1.5	13.0	10.0	44.0	Ø 34.0	706108
DN015	19.0 x 1.5	19.0	16.0	44.0	Ø 34.0	706109
DN020	23.0 x 1.5	23.0	20.0	44.0	Ø 50.5	706110
DN025	29.0 x 1.5	29.0	26.0	44.0	Ø 50.5	706111
DN032	35.0 x 1.5	35.0	32.0	48.0	Ø 50.5	706112
DN040	41.0 x 1.5	41.0	38.0	48.0	Ø 64.0	706113
DN050	53.0 x 1.5	53.0	50.0	49.0	Ø 77.5	706114
DN065	70.0 x 2.0	70.0	66.0	53.0	Ø 91.0	706115
DN080	85.0 x 2.0	85.0	81.0	57.0	Ø 106.0	706116
DN100	104.0 x 2.0	104.0	100.0	60.0	Ø 130.0	706117
DN125*	129.0 x 2.0	129.0	125.0	60.0	Ø 155.0	710629
DN150*	154.0 x 2.0	154.0	150.0	60.0	Ø 183.0	710630
DN200*	204.0 x 2.0	204.0	200.0	60.0	Ø 233.5	710631

ISO	Abmessung dimension	A	B	C	E	MB 1.4404
DN008	13.5 x 1.6	13.5	10.3	44.0	Ø 34.0	706546
DN010	17.2 x 1.6	17.2	14.0	44.0	Ø 34.0	706548
DN015	21.3 x 1.6	21.3	18.1	44.0	Ø 34.0	706550
DN020	26.9 x 1.6	26.9	23.7	44.0	Ø 50.5	706552
DN025	33.7 x 2.0	33.7	29.7	48.0	Ø 50.5	706554
DN032	42.4 x 2.0	42.4	38.4	48.0	Ø 64.0	706556
DN040	48.3 x 2.0	48.3	44.3	50.0	Ø 64.0	706558
DN050	60.3 x 2.0	60.3	56.3	54.0	Ø 91.0	706560
DN065	76.1 x 2.0	76.1	72.1	57.0	Ø 106.0	706562
DN080	88.9 x 2.3	88.9	84.3	62.0	Ø 119.0	706564
DN100*	114.3 x 2.3	114.3	109.7	62.0	Ø 144.4	706566

ASME BPE	Abmessung dimension	A	B	C	E	MB 1.4404
DN 1/2"	12.7 x 1.65	12.7	9.4	44.0	Ø 34.0	710632
DN 3/4"	19.05 x 1.65	19.05	15.75	44.0	Ø 34.0	710633
DN 1"	25.4 x 1.65	25.4	22.1	44.0	Ø 50.5	710634
DN 1 1/2"	38.1 x 1.65	38.1	34.8	48.0	Ø 64.0	710635
DN 2"	50.8 x 1.65	50.8	47.5	49.0	Ø 77.5	710636
DN 2 1/2"	63.5 x 1.65	63.5	60.2	55.0	Ø 91.0	710637
DN 3"	76.2 x 1.65	76.2	72.9	57.0	Ø 106.0	710638
DN 4"	101.6 x 2.11	101.6	97.38	59.0	Ø 130.0	710639

*) nicht in der DIN-Norm enthalten | non-standard