



5.6 MEDI-TK® DIN11852-TK

Hygiene T-Stück -kurz/egal-

hygienic tee -short/equal-

Legende

- A Außendurchmesser
- B Innendurchmesser
- C Gesamtlänge
- D Länge Abgang

Abmessungen

DIN 11866-A (DIN 11850)

Normative Verweisungen, Anforderungen / Abnahmen (Vormaterial)

DIN 11850
 DIN 11852
 EN 10088-1
 EN 10204

Kennzeichnung

Hersteller, Werkstoff, Abmessung, Schmelze

Toleranzen (Bezugsnorm)

DIN 11852 / ISO 2768-1

Verpackung

einzelnd verpackt

auf Wunsch:

Schweißenden mit Kunststoffkappen verschlossen

Werkstoffe

1.4404

weitere Werkstoffe auf Anfrage

Oberflächenbeschaffenheit ("MB/MB")

innen: metallblank, gebeizt
 (Vormaterial: Ra < 0.8/1.6 µm)
 außen: metallblank

Key

external diameter
 internal diameter
 total length
 length of the branch

Dimensions

DIN 11866-A (DIN 11850)

Normative References, Requirements / Approvals (raw material)

DIN 11850
 DIN 11852
 EN 10088-1
 EN 10204

Marking

manufacturer, material, dimension, heat

Tolerances (reference standard)

DIN 11852 / ISO 2768-1

Packaging

wrapped separately

on demand:

welding ends covered with plastic caps

Materials

316L

further materials on request

Surface condition ("MB/MB")

inside: bright finish, pickled
 (raw-material: Ra < 0.8/1.6 µm)
 outside: bright finish

MB/MB Innen- und Außenoberfläche unbearbeitet | raw internal and external surface finish

MB/K400 Innenoberfläche unbearbeitet, außen geschliffen K400 | raw internal surface, outside mech. polished K400




Oberflächenqualitäten | surface qualities

Innenoberfläche internal surface finish	(Zusatz-) Artikel-Nr. (additional) order No
e-polieret e-polished	999277
Ra < 0.8 µm (STERI-line®)	999284
Ra < 0.6 µm (STERI-line®)	999279
Ra < 0.5 µm (STERI-line®)	999280
Ra < 0.4 µm (HQ-line®)	999281
Ra < 0.25 µm (HQ-line®)	999283

Bestellbeispiel | ordering example

Artikelbezeichnung + Ausführungen article description + variants	Artikel-Nr. order No
MEDI-TK® DN025 DIN (1.4404) DIN 11852-TK, außen K400 hygienic tee -short/equal-, outside polished K400	711130
e-polieret e-polished	+ 999277
Ra < 0.6 µm	+ 999279



weitere Abmessungen auf Anfrage | further dimensions on request

DIN	Abmessung dimension	A	B	C	D		MB/MB 1.4404	MB/K400 1.4404		
DN010	13.0 x 1.5	XXX	XXX	XXX	XXX		XXX	XXX		
DN015	19.0 x 1.5	19.0	16.0	70.0	11.0		702835	711128		
DN020	23.0 x 1.5	23.0	20.0	80.0	13.0		702956	711129		
DN025	29.0 x 1.5	29.0	26.0	100.0	16.0		702996	711130		
DN032	35.0 x 1.5	35.0	32.0	110.0	19.5		703022	711131		
DN040	41.0 x 1.5	41.0	38.0	120.0	23.0		703067	711132		
DN050	53.0 x 1.5	53.0	50.0	140.0	30.0		702945	711133		
DN065	70.0 x 2.0	70.0	66.0	160.0	40.0		702955	711134		
DN080	85.0 x 2.0	85.0	81.0	180.0	47.5		704249	711135		
DN100	104.0 x 2.0	104.0	100.0	200.0	58.0		705060	711136		
DN125	129.0 x 2.0	129.0	125.0	375.0	74.0		705385	711137		
DN150	154.0 x 2.0	154.0	150.0	450.0	90.0		705386	711138		
DN200	204.0 x 2.0	204.0	200.0	600.0	115.0		711140	711139		

reduzierte T-Stücke sind ebenfalls erhältlich
reduced tees are also available

